## ACEWAND<sup>®</sup> compatible Electrosurgical Units (ESU)

The ACEWAND<sup>®</sup> device can be connected to a standard electrosurgical pencil.

The ACEWAND<sup>®</sup> device should be used with an electrosurgical unit utilizing a conductive adhesive patient return electrode.

The ACEWAND<sup>®</sup> device is compatible with an electrosurgical unit made by the following manufacturers:

Covidien, ValleyLab, Medtronic, ConMed, Aspen, Megadyne, Ethicon, Olympus, Aesculap, Braun, Bovie, Aaron, Birtcher, Ellman, Surgitron, Storz.